


Exposure Indicator Strips

the Andersen sterilization system



- ▶ **Exposure Indicator Strips**
Immediate assurance of gas exposure! Strips provide a strong color change, rapid confirmation of EO exposure to items at end of sterilization cycle. Convenient self-stick backing adheres to conventional paper or cloth wrapping, including Andersen Seal and Peel®. Use one strip, per package, per cycle. For use with Anprolene and EOGas systems only.

Exposure Indicator Strips AN85 200 units/box \$27.50

 **ANDERSEN PRODUCTS**
sterility experts

www.sterility.com
Andersen Products, Inc.
Health Science Park • 3202 Caroline Drive
Haw River, NC 27258-9564 USA
800-523-1276 • fax 336-376-8153

www.andcal.com
H. W. Andersen Products of California, Inc.
Health Science Park • 3151 Caroline Drive
Haw River, NC 27258-9789 USA
800-524-3455 • fax 336-376-3088

Distributed in Europe by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex CO15 4XA UK
1255-428-328

P/INT19 012003