


AN-86

Package Closure Indicator Tape



Package Closure Indicator Tape

Yellow stripes will turn blue immediately after exposure to Ethylene Oxide, providing proof of exposure. Use the blank area in the middle to keep track of date and lot information.



Call 1-800-524-3455



Actual size: 4 3/8" x 1/2"

Get visual assurance of gas exposure and seal your packages in one step! Stripes provide a strong color change, and confirmation of EtO exposure. Self-adhesive tapes adhere to conventional paper, CSR wrap, cloth wrapping or Andersen Seal and Peel®. For use with Anprolene® and EOGas™ systems only.

Package Closure Indicator Tape AN86
200 units/box \$27.50

ANPRO ANDERSEN
PRODUCTS
the sterility experts

www.sterility.com
Andersen Products, Inc.
Health Science Park • 3202 Caroline Drive
Haw River, NC 27258-9564 USA
800-523-1276 • fax 336-376-8153

www.andcal.com
H. W. Andersen Products of California, Inc.
Health Science Park • 3151 Caroline Drive
Haw River, NC 27258-9575 USA
800-524-3455 • fax 336-376-3088

Distributed in Europe by:
H. W. Andersen Products, Ltd.
Davy Road • Clacton-on-Sea
Essex CO15 4XA UK
1255-428-328

P/N649 082103